ROTA-ADENO MonlabTest®

MO-076010 20 TESTS

Immunochromatographic test for the simultaneous detection of Rotavirus and Adenovirus in feces

For *in vitro* use only. Store at 2 - 30°C.

PRINCIPLE

The Rota-Adeno chromatographic immunoassay is a procedure for the *in vitro* qualitative detection of Rotavirus and Adenovirus antigens in human stool. The test is based on the immunological capture of colored microbeads during its passage along a membrane on which the monoclonal antibody has been immobilized. The Rota-Adeno MonlabTest uses a combination of 1) monoclonal antibodies against VP6 antigen of group A rotavirus, conjugated to red latex particles and specific monoclonal antibodies to rotavirus on the membrane. 2) monoclonal antibodies against adenovirus hexon antigen (present in all adenovirus subtypes), conjugated to blue latex particles and specific monoclonal antibodies to adenovirus on the membrane.

In this test, the sample is previously treated with a sample diluent to extract the rotavirus and adenovirus antigens present in the stool. After the extraction, the extract is simply added to the reaction device.

When the extract flows through the test membrane, the colored particles migrate. In the event of a positive result, specific antibodies on the membrane will capture the colored particles.

Depending on the virus content of the sample, different colored lines will be visible. These lines are used for interpretation of results, following five-minute incubation at room temperature.

INTENDED USE

Rotaviruses are the main cause of acute gastroenteritis, particularly in children under two years of age. Its discovery in 1973 and its association with child gastroenteritis represented a breakthrough in the study of gastroenteritis not due to acute bacterial infection. Transmission occurs through the fecal-oral route; the incubation period is between 1 and 3 days.

Adenoviruses are the second leading cause of viral gastroenteritis in children (10 -15%); they can also cause respiratory diseases and depending on the serotype, diarrhea, conjunctivitis, cystitis, and others. At least 47 adenovirus serotypes have been identified and in all of them the hexon antigen is present. Serotypes 40 and 41 are associated with gastroenteritis. The main clinical symptom of gastroenteritis caused by adenovirus is diarrhea, for 9 to 12 days, also occurring with fever and vomiting.

MATERIALS PROVIDED

- 20 cassettes
- 20 dilution buffer (1,5mL)
- Instructions for use

MATERIALS NOT INCLUDED IN THE KIT

- Test tubes
- Adjustable pipettes
- Pipette tips

PRECAUTIONS

1. Patient samples (stool) may contain infectious agents and must be handled and disposed of as potentially hazardous biological material.
2. Buffer contains sodium azide as an antimicrobial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if signs of contamination or precipitation are present.
3. Do not eat, drink, smoke or prepare food in areas where reagents or samples are handled.
4. Wear disposable gloves when handling samples. Thoroughly wash the hands after finishing working.
5. Do not exchange components from one kit lot to another.
6. All reagents are for *in vitro* use exclusively.
7. Before use, allow kit components and samples to reach room temperature, as cold reagents and/or samples can reduce test functionality. Twenty to 30 minutes are recommended to reach room temperature.
8. Do not use kit components beyond their expiration date.
9. In case of package damage, the product may be used if none of its components has been damaged.
10. It is important to add the accurate sample amount (4-5 drops of suspension). If the amount is less than the required, the chromatography may not be completed as the sample doesn't reach the reaction area, if it is more, nonspecific lines may appear instead of red and blue lines.
11. Used product should be disposed of according to applicable laws.
12. Do not use the test if any colored lines appear in the result area before using the test.
13. It is very important to take the adequate stool amount: 30-50 mg of solid feces or $\approx 100 \mu\text{l}$ of liquid feces. **A sample excess prevents an adequate chromatography.**

STORAGE

The product can be stored at any temperature between 2 and 30°C. Its expiry date is printed on the envelope.

SAMPLES

Stool samples should be collected as soon as symptoms arise, before 8 days have elapsed. From this moment, the number of viral particles in stool decreases and they are less reactive.

Do not use samples collected in transport media or added with preservatives (such as formalin, SAF, PVA or similar) or enrichment media as their presence could interfere with the proper performance of the test.

The use of fresh untreated samples is highly recommended. Samples may be stored refrigerated (4°C approximately) for 1 or 2 days before their analysis. For a prolonged storage, samples must be frozen at -20°C without any previous preparation. If this occurs, samples should be completely thawed, brought to room temperature, and homogenized before the assay.

Pay special attention when analyzing hemorrhagic samples as they often give false positive results when the blood content is high. An indicator of this destabilization of the test is usually the alteration of the color in the bands.

PROCEDURE

1. Important: Take sample from three different sample sites at least, in order to get a sample as much representative as possible.
2. Unscrew the cap from the vial with caution in order not to spill the extraction buffer. With the end of the applicator **collect** a sufficient amount of stool (**30 - 50 mg**). If feces are liquid use a pipette to transfer 100 microliters into the vial.
3. Insert the applicator with the sample into the vial. Screw the cap well and **shake** vigorously to ensure a homogeneous mixture.
4. Remove the reaction device from the aluminium pouch. Discard the desiccant bag as it only serves to preserve the test from humidity and does not take part in test performance.
5. Break the top of the vial.
6. Add **4 or 5 drops** in the sample area of the reaction device (round window marked with an arrow). Do not add particulate matter with the liquid.
7. Wait **10 minutes**, read and interpret the results.

For samples with CARY BLAIR media:

- From DELTALAB (DeltaSwab Cary Blair): Add 1.5 mL of sample resuspended in Cary Blair medium to the vial with the sample dilution buffer included in the kit. Shake and continue with the procedure from point 5.
- From other sources (e.g. Thermofisher): Add 0.5 mL of sample resuspended in Cary Blair medium to the vial with the sample dilution buffer included in the kit. Shake and continue with the procedure from point 5

READING OF RESULTS

NEGATIVE: A single **GREEN** line appears in the results area, aligned with the letter "C" (control) marked on the casing. This line should always appear.

POSITIVE:

Rotavirus: A **GREEN** line and a second **PINK/RED** line appear in the results area, aligned with the legend "T2" (test) marked on the casing. Color intensity may vary depending on the antigen concentration present in the sample.

Adenovirus: A **GREEN** line and a second **BLUE** line appear in the results area, aligned with the legend "T1" (test) marked on the casing. Color intensity may vary depending on the antigen concentration present in the sample.


Rotavirus/Adenovirus: A **GREEN** line, a **RED** line, and a **BLUE** line appear in the results area, aligned with the legends "C", "T2" and "T1" (test) marked on the casing.

INVALID: If no green line appears, the test will be invalid due to inadequate procedure, reagent damage, or addition of improper sample amount. Repeat the test using a new reaction device.

Sometimes, in rare occasions, when testing different matrix samples, nonspecific lines have been detected, what indicates a destabilization of the test. It is a situation whose incidence has been considered when reporting test parameters. In this case the result has to be considered invalid, and it is advisable to repeat the assay with an alternative technique.

Any line that appears beyond 10 minutes because of the nature of the sample, is of no diagnostic value.

Interpretation code for the images on the strip.

NOTE: The indicative letters located in the plastic device work as a guide and help to the interpretation of the test. The test line may not exactly match with such mark.

Final diagnosis should not be based exclusively on a single test result. It should be based on the consistency of test results with other applicable data and the clinical symptoms.

QUALITY CONTROL

If no green line appears the test is invalid, due to either inadequate performance or to reagent deterioration. In this case, repeat the test strictly adhering to the working protocol outlined in this instruction sheet.

WARNING: The inclusion of a control with known results is recommended to ensure that the data obtained are correct.

LIMITATIONS OF THE PROCEDURE

- The test should only be used for the detection of Rotavirus and Adenovirus antigens in feces.
- The test is qualitative, and no quantitative interpretation of results should be made, relating to positive line intensity.
- Several hundred samples were assessed to ensure proper test performance. The correlation of results with other techniques (ELISA) was excellent. However, interferences in test performance cannot be ruled out.
- A sample excess may result in the appearance of lines of indefinite color, instead of red and blue lines, or the test and control lines may not appear due to an inadequate chromatography. The brown lines have no diagnostic value. In this event, the test must be repeated with a lower sample amount or by diluting the prepared extract.
- No cross-reaction with other viruses or substances has been observed during test evaluation. A negative result does not preclude a potential Rotavirus or Adenovirus infection. The significance of results must be assessed in connection with patient clinical symptoms.
- The analysis of some samples may lead to the occurrence of indefinite color nonspecific bands, caused in most cases by negative samples. When these lines of indefinite color appear, the test must be repeated. In the event the same result is obtained, the analysis by other analytical method is recommended. There are different reasons for this kind of problem; among them, the most common ones are the presence of blood in the stool sample, using samples from patients with anti-antibodies or, also too high sample concentration. Therefore, it is recommended to use samples free of blood and pay attention to these points.
- The test may lead to positive results in feces from patients to whom the oral solution RotaTeq vaccine has been administered until 15 days before the administration.

SENSITIVITY AND SPECIFICITY

From diverse assays with several observers in different conditions and with different lots, we have established a typical sensitivity of 31 ng/mL for rotavirus and 31 ng/mL for adenovirus, although amounts of 8 ng/mL for rotavirus and 4 ng/mL for adenovirus may be detected at 10 minutes.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

A product of the technique is taken as reference test ELISAs for the detection of each antigen, approved in the Japanese market. Therefore, the results of the evaluation are as follows:

		rota-Adeno MonlabTest	
		+	-
ELISA	+	94	1
	-	12	307

Concordance percentage: 100% [(94+307)/414] = 96.9%.

Positive concordance: 100% [94 / (94+1)] = 98.9%.

Negative concordance: 100% [307 / (307+12)] = 96.2%.

Sensitivity concordance estimation-PPV: 100% [94 / (94+12)] = 88.7%.

Specificity concordance estimation-NPV: 100% [307 / (307+1)] = 99.7%

A product of the technique is taken as reference test ELISAs for the detection of each antigen, approved in the Japanese market. Therefore, the results of the evaluation are as follows:

		rota-Adeno MonlabTest	
		+	-
ELISA	+	47	1
	-	6	360

Concordance percentage: 100% [(47+360)/414] = 98.3%.

Positive concordance: 100% [47 / (47+1)] = 97.9%.

Negative concordance: 100% [360 / (360+6)] = 96.2%.

Sensitivity concordance estimation-PPV: 100% [47 / (47+6)] = 88.7%.

Specificity concordance estimation-NPV: 100% [360 / (360+1)] = 99.7%

REPEATABILITY

Three duplicates of each concentration from the sensitivity curve are tested using three different lots and the same results are obtained.

REPRODUCIBILITY

INTER-DAY PRECISION: Using 1 lot of the product, ten duplicates of the sensitivity curve are performed throughout ten consecutive days. Only a difference of one ½ dilution is observed, acceptable and tolerable for the assay.

INTER-LAB PRECISION: Three different laboratories-operators test those same samples, presenting high precision and concordance. Only a difference of one ½ dilution is observed, acceptable and tolerable for the assay.

INTER-LOT PRECISION: Using 3 different lots of the product, a sensitivity curve is constructed in duplicate. The analysis is performed by the same person in the same day. Only a difference of one 1/2 dilution is observed, acceptable and tolerable for the assay.

HOOK EFFECT

No hook effect has been observed at a concentration equivalent to 65 times (2000 ng/mL) the limit of the test (31 ng/mL Rota + 31 ng/mL Adeno).

INTERFERING SUBSTANCES

The substances listed in the table, at the indicated concentration, didn't interfere with the results. Three lots were used for the study.

Substance	Initial concentration	Substance	Initial concentration
Bilirubin F	0,9 – 1,1 mg/mL	Atropine	(40 mg/dL)
Bilirubin C	0,9 – 1,1 mg/mL	Caffeine	(40 ng/dL)
Hemoglobin	22 – 27 mg/mL	Gentisic acid	(40 mg/dL)
Acetaminophen	(20 mg/dL)	Urea	(4.000 mg/dL)
Acetylsalicylic acid	(20 mg/dL)	Uric acid	(10 mg/dL)
Ampicillin	(40 mg/dL)	Ascorbic acid	(100 mg/dL)
Glucose	(2.000 mg/dL)		

INTERFERING MICROORGANISMS

The microorganisms listed, at the indicated concentration, did not interfere with the results.

Bacteria

Concentration (McF: 5(1x10⁸/mL)

<i>E.coli</i>	Groupe G Strep	<i>H.influenzae</i>
<i>K.pneumoniae</i>	<i>S.agalactiae</i>	<i>C.albicans</i>
<i>E.faecalis</i>	<i>S.mutans</i>	<i>B.catarrhalis</i>
<i>S.aureus</i>	<i>S.pneumoniae</i>	<i>B.pertussis</i>
<i>S.epidermidis</i>	<i>S.marcescens</i>	<i>P.vulgaris</i>
<i>S.haemolyticus</i>	<i>P.aeruginosa</i>	<i>S.saprophyticus</i>
Groupe B Strep	<i>N.gonorrhoeae</i>	<i>S.millieri</i>
Groupe C Strep		<i>E.faecium</i>

Virus-1

Virus	Subtype	Final concentration TCID ₅₀ /mL
Influenza type A	H1N1	1.4x10 ⁶
Influenza type A	H3N2	3.0x10 ⁶
Influenza type B	Victoria	1.6x10 ⁶
Influenza type B	Yamagata	1.5x10 ⁷

Virus-2

Virus	Subtype	Final concentration TCID ₅₀ /mL
Echovirus	6	7.2x10 ⁵
Poliovirus	2	2.2x10 ⁷
Parainfluenza	2	2.2x10 ⁷

REFERENCES

- F. Bon et al. *Prevalence of a group A rotavirus, human calicivirus, astrovirus, and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France.* Journal of Clinical Microbiology, Sept. 1999, p. 3055-3058
- Bodo R. Eing et al. *Evaluation of two enzyme immunoassays for detection of human rotaviruses in fecal specimens,* Journal of Clinical Microbiology, Dec. 2001, p. 4532-4534
- Umesh D. Parashar et al. *Global illness and deaths caused by rotavirus disease in children,* Emerging Infectious Diseases, vol. 9, No.5, May 2003, p. 565-572

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

Manufactured by

For in vitro diagnostic use

Do not re-use

Please read pack insert

Contains sufficient for <n> tests

Dry storage

Catalogue number

Store at

Lot number

Expiry date

This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices

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

- F. Bon et al. *Prevalence of a group A rotavirus, human calicivirus, astrovirus, and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France.* Journal of Clinical Microbiology, Sept. 1999, p. 3055-3058
- Bodo R. Eing et al. *Evaluation of two enzyme immunoassays for detection of human rotaviruses in fecal specimens,* Journal of Clinical Microbiology, Dec. 2001, p. 4532-4534
- Umeh D. Parashar et al. *Global illness and deaths caused by rotavirus disease in children,* Emerging Infectious Diseases, vol. 9, No.5, May 2003, p. 565-572

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