One-Step Amphetamine MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step Amphetamine MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of amphetamine in human urine at a cut-off of 1000ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

Amphetamines (amphetamine, methamphetamine, and the structurally related "designer" drugs, e.g. "Ecstasy") are sympathomimetic amines whose biological effects include potent central nervous system (CNS) stimulation, anorectic, hyperthermic, and cardiovascular properties. They are usually taken orally, intravenously, or by smoking. Amphetamines are readily absorbed from the gastrointestinal tract and are then either deactivated by the liver or excreted unchanged in the urine. Amphetamine is metabolized to deaminated (hippuric and benzoic acids) and hydroxylated metabolites. Methamphetamine is partially metabolized to amphetamine, its major active metabolite. Amphetamines increase the heart rate and blood pressure and suppress the appetite. Some studies indicate that heavy abuse may result in permanent damage to certain essential nerve structures in the brain. Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Amphetamine MonlabTest® is based on the principle of the highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of amphetamine in human urine at 1000ng/mL cut-off concentration.

PRINCIPLE

The One-Step Amphetamine MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which was pre-coated with drug conjugate on the test band.

A colored anti-amphetamine monoclonal antibody colloidal gold conjugate pad is placed at the right end of the membrane. In the absence of drug in the urine, the solution of colored antibody-colloidal gold conjugate and urine moves upward

chromatographically by capillary action across the membrane. This solution then migrates to the immobilized drug-protein conjugate zone to form a visible line as the antibody complexes with the drug-protein conjugate. Therefore, the formation of a visible precipitant in the test band region occurs when the test urine is **negative** for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the drug-protein conjugate on the test band region for the limited antibody sites on the anti-amphetamine antibody-colloidal gold conjugate. When a sufficient concentration of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein conjugate zone on the test band region (T). Therefore, absence of the color band on the test region indicates a **positive** result.

A control band that has a different antigen/antibody reaction is added to the immunochromatographic membrane strip at the control region (C) to indicate that the test has performed properly. This control line should always appear regardless of the presence of drug or metabolite. This means that negative urine will produce two colored bands, and positive urine will produce only one band. The presence of this colored band in the control region also serves as verification that 1) sufficient volume has been added, and 2) that proper flow was obtained.

MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each.
- One instruction sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer

STORAGE AND STABILITY

The test kit should be stored refrigerated or at temperature 2-30°C in the sealed pouch for the duration of the shelf-life.

PRECAUTIONS

- FOR IN VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.

SPECIMEN COLLECTION AND HANDLING

The One-Step Amphetamine MonlabTest® is formulated for use with urine specimens. Fresh urine does not require handling or pretreatment. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days, or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing.

Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.


Note: Urine specimens and all materials coming in contact with them should be handled and disposed of as if it is infectious and capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

TEST PROCEDURE


- Review "Specimen Collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient name or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15ml) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

INTERPRETATION OF RESULTS

Negative


Two colored lines should be observed in the viewing window. The line in the test region (T) is the test line; the line in the control region (C) is the control line, which is used to indicate proper performance of the device. The color intensity of the test line may be weaker or stronger than that of the control line.

Positive:

Only one colored line appears in the control region (C). The absence of a test line indicates a positive result.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new test device.

Note: A very faint line on the test region indicates that the amphetamine in the sample is near to the cut-off level for the test. These samples should be re-tested or confirmed with a more specific method before a positive determination is made.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well as other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** or lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Quality control specimens are available from commercial sources when testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of the One-Step Amphetamine MonlabTest® was evaluated in comparison to a commercially available immunoassay at a cut-off of 1000ng/mL for amphetamine. One hundred (100) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

In a separate study, seventy two (72) urine samples, obtained from a clinical laboratory where they were screened and confirmed positives by the commercially available immunoassay and GC/MS, were tested by One-Step Amphetamine MonlabTest®. The concentration of amphetamine in the urine samples ranged from 550ng/mL to 4600ng/mL. Of the fifty five (55) samples with amphetamine concentration ≥ 1500 ng/mL, all were found to be positive by both methods (100% agreement). Of the six (6) samples with amphetamine concentration from 1250 to 1500ng/mL, all were found to be positive by both methods (100% agreement). Of the eleven (11) samples with amphetamine concentration from 550 to 810ng/mL, all were found to be negative by both methods (100% agreement).

Reproducibility

The reproducibility of the One-Step Amphetamine MonlabTest® was evaluated at four different sites using blind controls. Of fifty (50) samples with amphetamine concentration of 500ng/mL, all were determined negative. Of the fifty (50) samples with amphetamine concentration of 2000ng/mL all were determined positive.

Precision

The precision of the One-Step Amphetamine MonlabTest® was determined by conducting the test with spiked controls. The control at the 500ng/mL should give a negative result and the control at the 1500ng/mL should give a positive result.

Concentration (ng/mL)	Number Tested	Correct Result	%Correct Result
500	50	50	100
1500	50	50	100

Specificity

The specificity for the One-Step Amphetamine MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
d-Amphetamine	1,000
l-Amphetamine	1,000
(+/-) 3, 4-methylenedioxyamphetamine (MDA)	5,000

The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.

Acetone	Acetaminophen
Albumin	Amitriptyline
Ampicillin	Aspartame
Aspirin	Atropine
Benzocaine	Bilirubin
Caffeine Chloroquine	(+)-Chlorpheniramine
(+/-)-Chlorpheniramine	Creatine pheniramine
Dextromethorphan Dexbrom	- Dimethylaminoantipyrine
Doamine	(-)-Epinephrine
Erythromycin	Ethanol
Furosemide	GuaiacolGlycerylEther
Glucose	Hemoglobin
Imipramine	(+/-)-Isoproterenol
(1R, 2S)-(-)-N-Methyl-Ephedrine	Lidocaine
(+/-)3,4-methylenedioxyamfetamine	(+)-Naproxen
etamine	
(+/-)-Norephedrine	Oxalic Acid
Penicillin-G	Pheniramine
Phenothiazine	l-Phenylephrine
Quinidine	Riboflavin
Sodium Chloride	Sulindac
Thioridazine	Trifluoperazine
Trimethobenzamide	Tyramine
Vitamin C	

REFERENCES

- Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn. M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts Illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

PACKAGING

CHEM-7003C-40

40 Amphetamine Tests

GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

Manufacturer: Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International Medical Park, Pudong, 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

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REFERENCIAS

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)