

## Haptoglobin MonlabTest®

Turbidimetry



### Quantitative determination of Haptoglobin (HAPTO)

**IVD** Only for professional in vitro diagnostic use  
Store at 2-8°C

#### INTENDED USE

HAPTO is a quantitative turbidimetric test for the measurement of haptoglobin in human serum or plasma.

#### PRINCIPLE OF THE METHOD

Anti-haptoglobin antibodies when mixed with samples containing haptoglobin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the haptoglobin concentration of the patient sample, that can be quantified by comparison from a calibrator of known haptoglobin concentration.

#### CLINICAL SIGNIFICANCE

The haptoglobin is an  $\alpha_2$ -glycoprotein synthesized in the liver that binds hemoglobin irreversibly. The hapto-hemoglobin complexes, as well as free haptoglobin itself, play significant roles in the iron storage and prevent possible renal damage as a consequence of hemoglobin excretion. As an acute-phase protein, haptoglobin is increased in the presence of acute inflammatory process, tissue necrosis or malignancy.

Haptoglobin deficiency in plasma is a consequence of hemolysis "in vivo", presence of estrogens in pregnancy and oral contraceptive therapy, as well as most forms of acute or chronic hepatocellular disease, including viral hepatitis.

Haptoglobin test is mainly used for the determination and monitoring of the hemolytic disorders. Under normal circumstances, approximately 1% of circulating red blood cells are destroyed every day. If this increases to 2%, it will completely deplete plasma haptoglobin in the absence of production stimulus such as acute inflammation or corticosteroids therapy.

#### REAGENTS

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Sodium azide 0.95 g/L.
<b>Antiserum (R2)</b>	Goat serum, anti-human haptoglobin pH 7.5. Sodium azide 0.95 g/L.
<b>Optional</b>	Multicalibrator Protein Serum (MO-165044).

#### CALIBRATION

The assay has been standardized against the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements, IRMM). It must be used the Multicalibrator Protein Serum Calibrator to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

#### PREPARATION

- Reagents: Ready to use.
- Calibration Curve: Prepare the following Multicalibrator Protein Serum dilutions in CIna 9 g/L as diluent. Multiply the concentration of the haptoglobin calibrator by the corresponding factor stated in table below to obtain the haptoglobin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator ( $\mu$ L)	--	10	25	50	75	100
NaCl 9 g/L ( $\mu$ L)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

#### ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 340 nm filter (320 – 360 nm).

#### SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.  
The samples with presence of fibrin should be centrifuged before testing.  
The samples with presence of fibrin should be centrifuged before testing.  
Do not use highly hemolyzed or lipemic samples.

#### PROCEDURE

1. Bring the reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions:  
Wavelength: 340 nm  
Temperature: 37 °C  
Cuvette light path: 1cm
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette:

Reagent R1 ( $\mu$ L)	800
Sample or Calibrator ( $\mu$ L)	10

5. Mix and read the absorbance ( $A_1$ ) after the sample addition.
6. Immediately, pipette into the cuvette:

Reagent R2 ( $\mu$ L)	200
-----------------------	-----

7. Mix and read the absorbance ( $A_2$ ) of calibrators and sample exactly 2 minutes after the R2 addition.

**MONLAB has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

#### CALCULATIONS

Calculate the absorbance difference ( $A_2 - A_1$ ) of each point of the calibration curve and plot the values obtained against the haptoglobin concentration of each calibrator dilution. Haptoglobin concentration in the sample is calculated by interpolation of its ( $A_2 - A_1$ ) in the calibration curve.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. MONLABTEST Multicontrol Protein Serum (MO-165045) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### REFERENCE VALUES<sup>2</sup>

Between 30 - 200 mg/dL. Each laboratory should establish its own reference range.

#### PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 300 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 1.3 mg/dL give non-reproducible results.



