Calibrador HDL_c/LDL_c
MonlabTest®
 Suero Humano

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HDL_c/LDL_c Calibrator
MonlabTest®
 Human Serum

Quantitative determination of HDL_c/LDL_c

Only for professional *in vitro* diagnostic use.
 Store at 2-8°C.

PRODUCT CHARACTERISTICS

HDL_c/LDL_c calibrator is a human lyophilised serum.
 Is intended for determination of HDL_c or/and LDL_c concentrations with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives.
 The concentration of the components is lot-specific.
 See chart below for the values of ingredients.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Information available on request.

PREPARATION

-Open one vial carefully, without loss of the lyophilised.
 -Reconstitute (→) the lyophilised with 1,0 mL of distilled water.
 -Mix thoroughly, avoiding foam forming.
 -Bring to room temperature for about 30 minutes before use.
 Improper handling and/or storage can affect results.
 Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial grown.

Store tightly capped when not use.

- After reconstitution is stable for:

- At 15°C to 25°C 30 hours
- At 2°C to 8°C 14 days
- At -25°C to -15°C 3 months

PRESENTATION

MO-165111

4 x 1 mL

LOT

Component	Method	Value	Units
HDL _c cholesterol	Direct- Liquid		mg/dL mmol/L
LDL _c cholesterol	Direct- Liquid		mg/dL mmol/L