
<p><b>FR Calibrador</b> <b>MonlabTest®</b> Turbidimetría</p>	 
--	--

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**Quantitative determination of Rheumatoid Factors (RF)**

Only for professional *in vitro* diagnostic use  
Store at 2-8°C

**PRODUCT CHARACTERISTICS**

The calibrator is a lyophilised serum from human origin containing a RF concentration suitable to calibrate RF determinations by turbidimetry method.

**REAGENTS**

<b>RF Calibrator</b>	Calibrator. Human serum. RF concentration is stated on the vial label.
----------------------	--

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

The concentration value of the Calibrator has been standardized against the International Reference Standard from NIBSC 64/002.

**PREPARATION**

Reconstitute (→) with 2,0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring to room temperature for about 10 minutes before use.

**STORAGE AND STABILITY**

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8°C or 3 months at -20°C.

**PROCEDURE**

To be used in turbidimetric assays.

**PRESENTATION**

MO-165047	1 x 2 mL
-----------	----------