## β<sub>2</sub>-Microglobulina MonlabTest®

Turbidimetría Látex.

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1. Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
2. Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
3. Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
4. Wibell L et al. Nephron 1973; 10: 320-331.
5. Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
6. Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
7. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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IVD

### Quantitative determination of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m)

Only for professional *in vitro* diagnostic use.

Store at 2-8°C.

#### PRINCIPLE OF THE METHOD

The β<sub>2</sub>-microglobulin MonlabTest is a quantitative turbidimetric test for the measurement of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m) in human serum or plasma. Latex particles coated with anti-human β<sub>2</sub>-m are agglutinated when mixed with samples containing β<sub>2</sub>-m. The agglutination causes an absorbance change, dependent upon the β<sub>2</sub>-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

#### CLINICAL SIGNIFICANCE

The β<sub>2</sub>-microglobulin is a protein located on the surface of human lymphocytes and other nucleated cells. Free β<sub>2</sub>-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β<sub>2</sub>-m is a sensitive indicator of renal insufficiency. Also, the β<sub>2</sub>-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

#### REAGENTS

<b>β<sub>2</sub>-m Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8.2. Preservative.
<b>β<sub>2</sub>-m Latex (R2)</b>	Particles coated with goat IgG anti-human β <sub>2</sub> -m, pH 7.5. Preservative.
<b>β<sub>2</sub>-m CAL</b>	Calibrator. β <sub>2</sub> -m concentration is stated on the vial label.
<b>Optional:</b>	β <sub>2</sub> -m Control MonlabTest (MO-165054).

#### PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product. Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

#### CALIBRATION

Use β<sub>2</sub>-Microglobulin Calibrator MonlabTest (MO-165053). The sensitivity of the assay and the target value of the calibrator have been standardized against the 1<sup>st</sup> International ERM-DA470k/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

#### PREPARATION

##### β<sub>2</sub>-m Calibrator:

Serum method: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

**Calibration curve:** Prepare the following β<sub>2</sub>-m calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the β<sub>2</sub>-m calibrator by the corresponding factor stated in table below to obtain the β<sub>2</sub>-m concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
β <sub>2</sub> -m Calibrator (μL)	--	5	10	25	50	100
NaCl 9 g/L (μL)	100	95	90	75	50	--
Factor	0	0,05	0,1	0,25	0,5	1,0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

**Do not freeze.** Frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

**β<sub>2</sub>-m Calibrator:** Stable for 1 month at 2-8°C or 3 months at -20°C.

#### ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostable at 37°C with a 540 nm filter.

#### SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with particles or fibrin should be centrifuged before testing. Do not use hemolyzed or lipemic samples.

#### PROCEDURE

1. Bring the Reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions:  
Wavelength: ..... 540 nm (530-550)  
Temperature: ..... 37°C  
Cuvette light path: ..... 1 cm
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette:

Diluent R1	800 μL
Latex R2	200 μL
Calibrator or sample	10 μL

5. Mix and read the absorbance immediately (A<sub>1</sub>) and after 3 minutes (A<sub>2</sub>) of the sample addition.

**MONLAB has instruction sheets available for several automatic analyzers. Instructions for many of them are available on request.**

#### CALCULATIONS

Calculate the absorbance difference (A<sub>2</sub>-A<sub>1</sub>) of each point of the calibration curve and plot the values obtained against the β<sub>2</sub>-m concentration of each calibrator dilution. β<sub>2</sub>-m concentration in the sample is calculated by interpolation of its (A<sub>2</sub>-A<sub>1</sub>) in the calibration curve.

#### QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the β<sub>2</sub>-m Control MonlabTest (MO-165054). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### REFERENCE VALUES

Serum: from 1.0 to 3.0 mg/L.  
Each laboratory should establish its own reference range.

#### PERFORMANCE CHARACTERISTICS

1. **Linearity limit:** from quantitation limit of 0.150 mg/L to the upper calibration level. If the results obtained are greater than the linearity limit, then dilute the sample 1/5 with NaCl 9 g/L and retest again.
2. **Detection limit:** values less than 0.100 mg/L give non-reproducible results.
3. **Prozone effect:** no prozone effect was detected up to 100 mg/L.
4. **Sensitivity:** Δ 0.048 A. mg/L.
5. **Precision:** The reagent has been tested for 20 days, using three different β<sub>2</sub>-m concentrations in an EP5-based study.

EP5	CV (%)		
	+/- 1 mg/L	+/- 3.2 mg/L	+/- 8.5 mg/L
Total	4.0%	3.4%	1.7%
Within Run	2.8%	2.0%	1.2%
Between Run	1.7%	1.5%	1.2%
Between Day	2.2%	2.4%	0.0%

6. **Accuracy:** results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. Forty-two samples of different β<sub>2</sub>-m concentrations were assayed. The correlation coefficient (r) and the regression equation were:  
y = -0.0637 + 0.986 x; r = 0.999.  
The results of the performance characteristics depend on the analyzer used.

#### INTERFERENCES

Bilirubin (20 mg/L), hemoglobin (10 g/L) and lipids (10 g/L) do not interfere. Rheumatoid factors (150 IU/mL) interfere. Other substances may interfere<sup>7</sup>.

#### BIBLIOGRAPHY

1. Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
2. Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
3. Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
4. Wibell L et al. Nephron 1973; 10: 320-331.
5. Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
6. Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
7. Young DS. Effects of drugs on clinical laboratory test, 4<sup>th</sup> ed. AACC Pres, 1995.

#### PACKAGING

MO-165031  
R1: 1 x 40 mL  
R2: 1 x 10 mL  
β<sub>2</sub>-m Calibrator: 1x1 mL

#### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by