**HbA<sub>1c</sub> MonlabTest®**

Turbidimetría Látex.

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**HbA<sub>1c</sub> MonlabTest®**  
Latex Turbidimetry

**Quantitative determination of glycosylated hemoglobin (HbA<sub>1c</sub>) in human blood.**  
Only for professional *in vitro* diagnostic use. Store at 2-8°C.

**PRINCIPLE OF THE METHOD**

This method utilizes the interaction of antigen and antibody to directly determine the HbA<sub>1c</sub> in whole blood. Total hemoglobin and HbA<sub>1c</sub> have the same unspecific absorption rate to latex particles. When mouse antihuman HbA<sub>1c</sub> monoclonal antibody is added (R2), latex-HbA<sub>1c</sub>-mouse anti human HbA<sub>1c</sub> antibody complex is formed. Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA<sub>1c</sub> absorbed on to the surface of latex particles. The amount of agglutination is measured as absorbance. The HbA<sub>1c</sub> value is obtained from a calibration curve.

**CLINICAL SIGNIFICANCE**

Throughout the circulatory life of the red cell, Hemoglobin A<sub>1c</sub> is formed continuously by the addition of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al<sup>1</sup> showed Hemoglobin A<sub>1c</sub> in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that Hemoglobin A<sub>1c</sub> serve as an indicator of metabolic control of the diabetic, since Hemoglobin A<sub>1c</sub> levels approach normal values for diabetics in metabolic control.<sup>2,3,4</sup> Hemoglobin A<sub>1c</sub> has been defined operationally as the "fast fraction" hemoglobins (HbA<sub>1a</sub>, A<sub>1b</sub>, A<sub>1c</sub>) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA<sub>0</sub>. The present procedure utilizes an antigen and antibody reaction to directly determine the concentration of the HbA<sub>1c</sub>.

**REAGENTS**

<b>R1</b>	Latex 0.13%, Buffer, stabilizer.
<b>R2</b>	Mouse anti-human HbA <sub>1c</sub> monoclonal antibody 0.05mg/mL, goat anti-mouse IgG polyclonal antibody 0.08mg/dL, Buffer, stabilizers.
<b>R3 (Hemolysis reagent)</b>	Water and stabilizers
<b>Optional</b>	HbA <sub>1c</sub> Calibrator (4 levels) MonlabTest – (MO-165061) HbA <sub>1c</sub> Control (2 levels) MonlabTest – (MO-165062)

**PRECAUTIONS**

All human specimens should be regarded as potentially biohazardous. Therefore, universal precautions should be used in specimen handling (gloves, lab garments, avoid aerosol production, etc.).

**PREPARATION**

R1, R2 and R3 are ready to use. Mix gently before use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

R1 and R2 are stable for at least one month after opening stored at 2-8°C.

**Reagent deterioration:** Alterations in the physical appearance of the reagents or values of control materials outside of the manufacturer's acceptable range may be an indication of reagent instability.

**ADDITIONAL EQUIPMENT**

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 660 nm filter.
- General laboratory equipment (Note 1)

**SAMPLES**

Special preparation of the patient is unnecessary. Fasting specimens are not required. No special additives or preservatives other than anticoagulants are required. Collect venous blood with EDTA using aseptic technique. HbA<sub>1c</sub> in whole blood collected with EDTA is stable for one week at 2-8°C.<sup>5</sup>

To determine HbA<sub>1c</sub>, a hemolysate must be prepared for each sample:

1. Dispense 1 mL Hemolysis Reagent into tubes labeled: Calibrator, Control, Patients, etc. Note: Plastic or glass tubes of appropriate size are acceptable.
2. Place 20 µL of well mixed whole blood into the appropriately labeled lyse reagent tube. Mix.
3. Allow to stand for 5 minutes or until complete lysis is evident. Hemolysates may be stored up to 10 days at 2-8°C

**PROCEDURE**

1. Assay conditions:  
Wavelength: ..... 660 nm (600-660)  
Temperature: ..... 37°C  
Cuvette lighth path: ..... 1 cm
2. Adjust the instrument to zero with distilled water.
3. Pipette into a cuvette (Note 2):

R1 (µL)	360
Calibrator or sample (µL)	10

4. Mix and incubate 5 minutes.
5. Pipette into the cuvette:

R2 (µL)	120
---------	-----

6. Mix and read the absorbance after 5 minutes (A) of the R2 addition.

**CALCULATIONS**

**HbA<sub>1c</sub> concentration (%)**

Plot (A) obtained against the HbA<sub>1c</sub> concentration of each calibrator (1 to 4 Level). HbA<sub>1c</sub> percentage in the sample is calculated by interpolation of its absorbance (A) in the calibration curve.

**QUALITY CONTROL**

HbA<sub>1c</sub> Control MonlabTest (MO-165062) is recommended to monitor the performance of manual and automated assay procedures. **Controls require hemolysis pretreatment after being reconstituted.**

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES<sup>11</sup>**

Recommended Values: less than 6% for a non-diabetic, less than 7% for glycemic control of a person with diabetes.

Each laboratory should establish its own expected values. In using Hemoglobin A<sub>1c</sub> to monitor diabetic patients, results should be interpreted individually. That is, the patient should be monitored against him or herself. There is a 3-4 week time lag before Hemoglobin A<sub>1c</sub> reflects changes in blood glucose level.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 2% to linearity limit of 16%.%

**Precision:**

Mean (%)	Intra-assay (n=20)		Inter-assay (n=20)	
	5,95	12,15	5,97	12,21
SD	0,19	0,18	0,14	0,15
CV (%)	3,20	1,47	2,31	1,24

**Sensitivity:** 1% = 0,056 (A)

**Accuracy:** Results obtained using MonlabTest reagents (y) did not show differences when compared with another commercial reagent (x). The results obtained using 40 samples were the following:

Correlation coefficient (r)<sup>2</sup> 0,995

Regression equation: y = 0,989x - 0,047

The results of the performance characteristics depend on the analyzer used.

**INTERFERENCES AND LIMITATIONS**

1. Bilirubin to 50mg/dL, ascorbic acid to 50mg/dL, triglycerides to 2000mg/dL, carbamylated Hb to 7.5mmol/L and acetylated Hb to 5.0mmol/L do not interfere in this assay.
2. It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin.<sup>6,7,8,9</sup>
3. It has been reported that elevated levels of HbF may lead to underestimation of HbA<sub>1c</sub> and, that uremia does not interfere with HbA<sub>1c</sub> determination by immunoassay.<sup>10</sup> It has been reported that labile intermediates (Schiff base) are not detected and therefore, do not interfere with HbA<sub>1c</sub> determination by immunoassay.<sup>5</sup>
4. It has been determined that Hemoglobin variants HbA<sub>2</sub>, HbC and HbS do not interfere with this method.
5. Other very rare variants of hemoglobin (e.g. HbE) have not been assessed.

**NOTES**

1. In order to avoid contamination, it is recommended to use disposable material.
2. Use clean disposable pipette for its dispensation.
3. MONLAB instruction sheets for several automatic analyzers.

**BIBLIOGRAPHY**

1. Trivelli, L.A., Ranney, H.M., and Lai, H.T., New Eng. J. Med. 284,353 (1971).
2. Gonen, B., and Rubenstein, A.H., Diabetologia 15, 1 (1978).
3. Gabbay, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J. Clin. Endocrinol. Metab. 44, 859 (1977).
4. Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978).
5. Tietz, N.W., Textbook of Clinical Chemistry, Philadelphia, W.B. Saunders Company, p.794-795 (1999).
6. Ceriello, A., et al, Diabetologia 22, p. 379 (1982).
7. Little, R.R., et al, Clin. Chem. 32, pp. 358-360 (1986).
8. Fluckiger, R., et al, New Eng.J. Med. 304 pp. 823-827 (1981).
9. Nathan, D.M., et al, Clin. Chem. 29, pp. 466-469 (1983).
10. Engbaek, F., et al, Clin. Chem. 35, pp. 93-97 (1989).
11. American Diabetes Association: Clinical Practice Recommendations (Position Statement). Diabetes Care 24 (Suppl. 1): S33-S55, (2001).

**PACKAGING**

MO-165035	R1: 1 x 30 mL
	R2: 1 x 10 mL
	R3: 1 x 125 mL

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS**

- |  |                                   |  |                                         |
|--|-----------------------------------|--|-----------------------------------------|
|  | Manufacturer                      |  | For <i>in vitro</i> diagnostic use only |
|  | Don't re-use                      |  | Consult instructions for use            |
|  | Contains sufficient for <n> tests |  | Keep dry                                |
|  | Catalogue Code                    |  | Temperature limitation                  |
|  | Lot Number                        |  | Use by                                  |