## One-Step Benzodiazepine MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

### INTENDED USE

The One-Step Benzodiazepines (BZO) MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of benzodiazepines in human urine at a cut-off of 300ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

### SUMMARY

Benzodiazepines are the most widely used anxiolytic drugs. They are used extensively as anti-anxiety agents, hypnotics, muscle relaxants and anticonvulsants. They are taken orally or sometimes by injection. Benzodiazepines also exhibit pharmacological activities. Benzodiazepines and metabolites are excreted in the urine. Their use can result in drowsiness, confusion. Benzodiazepines potentiate alcohol and other CNS depressants. Psychological and physical dependence on benzodiazepines can develop if high doses of the drug are given over a prolonged period.

Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Benzodiazepines MonlabTest® is based on the principle of the highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of benzodiazepines in human urine at a 300ng/mL cut-off concentration.

### PRINCIPLE

The One-Step Benzodiazepines (BZO) MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug-protein conjugate) competes with the drug which may be present in urine for limited antibody binding sites.

The test device contains a membrane strip which has been pre-coated with drug-protein conjugate on the test band region. A colored anti-benzodiazepines antibody-colloidal gold conjugate pad is placed at the end of the membrane. The colored antibody-colloidal gold conjugate moves along with urine, chromatographically by capillary action, across the membrane.

In the absence of drug in the urine, the colored antibody colloidal gold conjugate attaches to the drug-protein conjugate on the test band region to form a visible line as the antibody/drug-protein conjugate complexes. Therefore, the formation of a visible precipitant in the test band region occurs when the test urine is **negative** for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the limited antibody sites. When a sufficient amount of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the

colored antibody-colloidal gold conjugate to the drug-protein conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a **positive** result.

A control or reference band with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip to indicate that the test is performed properly. This control line should always appear regardless of the presence of drug or metabolite. This means that **negative** urine will produce **two** colored bands, and positive urine will produce only **one** band. The presence of this colored region also serves as verification that: 1) Sufficient volume has been added and 2) that proper flow was obtained.

### MATERIALS SUPPLIED

- 40 individually wrapped cassettes which include one disposable pipette each. Each test contains a membrane coated with drug-protein conjugate and a colloidal gold conjugate pad coated with anti-benzodiazepines antibody.
- One instruction sheet.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer

### STORAGE AND STABILITY

The test kit should be stored refrigerated or at temperature 2-30°C. Each device should remain in its sealed pouch for the duration of the shelf-life.

### PRECAUTIONS

- FOR IN-VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after expiration date.

### SPECIMEN COLLECTION AND HANDLING


The One-Step Benzodiazepines (BZO) MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

**Note:** Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

### TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15mL) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- **Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

## INTERPRETATION OF RESULTS


### Positive

The control line appears in the test window, but the test line is not visible.

### Negative

Two pink-rose lines (bands) are visible in the control (C) and test (T) areas of the test window. The intensity of the test line may be less than that of the control line; this still means negative result.

### Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

**Note:** A very faint line on the test region indicates that the benzodiazepines in the samples is near to the cut-off level for the test. These samples should be re-tested or confirmed with a more specific method before a positive determination is made.

## LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well as other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

## QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Before using a new kit with patient specimens, positive and negative controls should be tested. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with urine specimen.

## PERFORMANCE CHARACTERISTICS

### Accuracy

The accuracy of the One-Step BZO MonlabTest® was evaluated in comparison to a commercially available one-step immunoassay at a cut-off of 300ng/mL for benzodiazepines. One hundred (100) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

In a separate study, fifty three (53) urine samples, (obtained from a clinical laboratory where they were confirmed as positives by GC/MS) were tested with the One-Step Benzodiazepines MonlabTest® and a commercially available one-step immunoassay. Of the forty nine (49) samples with benzodiazepines concentration  $\geq 300$ ng/mL, all were found

to be positives by both methods (100% agreement). Of the four (4) samples with benzodiazepines concentration from 100 to 200ng/mL, all were determined negatives by both methods (100% agreement).

### Reproducibility

The reproducibility of One-Step BZO MonlabTest® was evaluated at four different sites using blind controls. Of the sixty (60) samples without oxazepam all were determined negative. Of the sixty (60) samples with oxazepam concentration of 600ng/mL, all were determined positive.

### Precision

The precision of the One-Step BZO MonlabTest® was determined by conducting the test with spiked controls. The control at 600ng/mL should give a positive result.

Concentration (ng/mL)	Number Tested	Correct Result	Correct Result %
0	50	50	100
600	50	50	100

### Specificity

The specificity for the One-Step BZO MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
Alprazolam	150
Bromazepam	800
Chlordiazepoxide	300
Clonazepam	1,000
Clobazam	200
Delorazepam	300
Clorazepam	100
Flunitrazepam	1,000
Diazepam	150
Lorazepam	1,500
Estazepam	150
Lormetazepam	1,000
Estazolam	2,500
Flurazepam	300
Medazepam	2,000
Nordiazepam	100
Nitrazepam	1,000
Oxazepam	300
Prazepam	1,000
Temazepam	150
Triazolam	1,500

**The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.**

Acetaminophen	Doxylamine
Acetone	(+/-)-Ephedrine
Albumin	(-)-Ephedrine
Amitriptyline	(+)-Epinephrine
Ampicillin	Erythromycin
Aspartame	Ethanol
Aspirin	Furosemide
Atropine	Glucose
Benzocaine	Guaicol glyceryl ether
Benzoylgonine	Hemoglobin
Bilirubin	Hydromorphone
Caffeine	hydroxy tyramine
(+)-Chlorpheniramine	(+/-)-Chlorpheniramine

Imipramine	Chlorpromazine
(+/-)-Isoproterenol	Creatine
Lidocaine	
Deoxyephedrine	Naloxone
Dexbrompheniramine	Naltrexone
Dextromethorphan	Naphthalene acetic acid
4-Dimethylaminoantipyrine	(+)-Naproxen
Dopamine	(+/-)-Norephedrine
11-not- $\Delta^9$ -THC-9-carboxylic acid	Quinidine
Oxalic Acid	Riboflavin
Oxycodone	Secobarbital
Penicillin-G	Sodium Chloride
Pentamine	Sulindac
Pentobarbital	Tenocyclidine
Pheniramine	Thiordazine
Phenobarbita	Trimethobenzamidel
Phenothiazine	Trifluoperazine
L-Phenylephrine	Tyramine
d-Propoxyhene	Vitamin C

**Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

**European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

Tel: 0049-40-2513175 / Fax: 0049-40-255726

Importer /

Distributor:

MONLAB SL

Cobalto, 74

08940 Cornellà de Llobregat (Spain)

Tel +34 934 335 860 / Fax +34 934 363 894

### REFERENCES

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Greenblatt, D, J, Shader, R,1, Benzodiazepines in Clinical Practice. New York Raven Press, 1974.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

### PACKAGING

CHEM-7010C-40

40 Benzodiazepines Cassettes

### GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

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## REFERENCIAS

- Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn. M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR (ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Greenblatt, D, J, Shader, R,1, Benzodiazepina in Clinical Practice. New York Raven Press, 1974.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.


- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-4OB, 1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.