Astrovirus MonlabTest®

MO-804009 20 TESTS

One step test to detect Astrovirus antigens

A rapid, one step test for the qualitative detection of Astrovirus antigens in stool samples.
For professional *in vitro* diagnostic use only.

INTENDED USE

The Astrovirus MonlabTest® is a rapid chromatographic immunoassay for the qualitative detection of Astrovirus antigens in human feces specimens to aid in the diagnosis of Astrovirus infection.

SYNTHESIS

Viral gastroenteritis is an infection caused by a variety of viruses that results in vomiting or diarrhea. Many different viruses can cause gastroenteritis, including rotaviruses, noroviruses, adenoviruses, sapoviruses, and astroviruses.

The main symptoms of viral gastroenteritis are watery diarrhea and vomiting. The affected person may also have headache, fever, and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness. Some research studies have shown that the duration of the symptoms are approximately three to four days. Astrovirus infection is not usually a severe situation and only in some rare cases leads to dehydration. Adenoviruses and astroviruses cause diarrhea mostly in young children, but older children and adults can also be affected.

PRINCIPLE

The Astrovirus MonlabTest® is a qualitative immunoassay for the detection of Astrovirus antigen in human feces samples. The membrane is pre-coated with monoclonal antibodies against Astrovirus antigens on the test line region. During testing, the sample reacts with the particle coated with anti-Astrovirus antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugate and generate a coloured line. A green coloured band always appears in the control line and serves as verification that sufficient volume was added, that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use the test if pouch is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test should be discarded in a proper biohazard container after testing.
- The test must be carried out within 2 hours of opening the sealed bag.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS PROVIDED

- 20 Test
- Instruction for use
- 20 specimen collection vial with buffer

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

SPECIMEN COLLECTION AND PREPARATION

Collect sufficient quantity of feces (1-2 g or mL for liquid sample). Stool samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-4°C/36-40°F) for 1-2 days prior to testing. For longer storage the specimen must be kept frozen at -20°C/-4°F. In this case, the sample will be totally thawed, and brought to room temperature before testing.

PROCEDURES

To process the collected stool samples (see illustration 1):

Use a separate specimen collection vial for each sample (with 1mL of buffer). Introduce the stick two times into the fecal specimen to pick up quite a lot of sample (125mg). Close the vial with the buffer and stool sample. Shake the vial in order to assure good sample dispersion. For liquid stool samples, aspirate the fecal specimen with a dropper and add 125µL into the specimen collection vial with buffer.

Test Procedure (see illustration 2)

Allow the tests, stool samples and buffer to reach to room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the Astrovirus MonlabTest® from its sealed pouch and use it as soon as possible.
2. Shake the specimen collection vial to assure good sample dispersion. Break off the cap of the vial.
3. Use a separate device for each sample. Dispense 4 drops into the specimen well (S). Start the timer.
4. Read the result at **10 minutes** after dispensing the sample.

Illustration 1




Illustration 2

INTERPRETATION OF RESULTS

Illustration 3

POSITIVE: Two lines appears across the central window in the result line region (**red** test line marked in the illustration 3 with the letter T) and in the control line region (**green** control line marked in the illustration 3 with the letter C).

NEGATIVE: Only one **green** band appears across the control line region marked in the illustration 3 with the letter C (control line).

INVALID: A total absence of the **green** control coloured band regardless the appearance or not of the **red** test line. Note: Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test:

- A green line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. Astrovirus MonlabTest® will only indicate the presence of Astrovirus in the specimen (qualitative detection) and should be used for the detection of Astrovirus antigens in feces specimens only. Neither the quantitative value nor the rate of increase in Astrovirus antigens concentration can be determined by this test.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of Astrovirus infection.
4. After one week of infection, the number of virus in feces decreases, making the sample less reactive. The stool samples should be collected during the first week of the presence of symptoms appears.
5. This test provides a presumptive diagnosis of Astrovirus infections. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

EXPECTED VALUES

A study of intestinal disease in the UK, published in 1999 determined incidence as 3.8/1000 patient years in the community (95%CI, range 2.3-6.4), the fourth most common known cause of viral gastroenteritis. Studies in the USA have detected astroviruses in the stools of 2-9% of children presenting symptoms; illness is most frequent in children of less than two years, although outbreaks among adults and the elderly have been reported. Early studies carried out in Glasgow demonstrated that a significant proportion of babies excreting virus particles, 12%, did not exhibit gastrointestinal symptoms, and seroprevalence studies carried out in the US have shown that 90% of children have antibody to HastV-1 by age 9, suggesting that (largely asymptomatic) infection is common. There is, as with most viral causes of gastroenteritis, a peak of incidence in the winter.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The evaluation was conducted comparing the results obtained using the Astrovirus MonlabTest® to a commercial available Astrovirus ELISA assay.

The results showed >94% of sensitivity and >99% of specificity.

Cross-Reactivity



It was performed an evaluation to determine the cross reactivity of Astrovirus MonlabTest®. There is not cross reactivity with common gastrointestinal pathogens, other organisms and substances occasionally present in feces.

- Rotavirus
- Adenovirus
- *Escherichia coli*
- *Campylobacter*
- *Giardia lamblia*
- Human Hemoglobin

REFERENCES

1. SUNITA SHASTRI et al., "Prevalence of Astroviruses in a Children's Hospital", JOURNAL OF CLINICAL MICROBIOLOGY Sept. 1998, p. 2571-2574 Vol. 36, No. 9.
2. ASHLEY et al., "Astrovirus-associated gastroenteritis in children", Journal of Clinical Pathology., 1978, 31, 939-943.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

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